

Proximate analysis of Coal

Introduction

Proximate Analysis includes analysis of basic things present in coal sample. It includes -

- Determination of % Moisture
- Determination of % Ash
- Determination of % Volatile matter
- Determination of % Fixed carbon

a. Determination of % Moisture

Procedure:

1. Weigh a clean dry and empty crucible (W_1 gm)
2. Take about 1gm of coal sample in the same crucible and weigh it again (W_2 gm).
3. Then place the silica crucible in muffle furnace or air oven, maintained at temperature range of $105-110^\circ\text{C}$ for 1 hour.
4. Cool it in a desiccator and weigh it again (W_3 gm). The loss in weight corresponds to moisture.

Observations:

Observations	Value
Weight of empty crucible (W_1)	= gm
Weight of crucible with coal sample (W_2)	= gm
Weight of crucible with coal sample after heating (W_3)	= gm
Weight of coal sample	$(W_2 - W_1) =$ gm

Formula and Calculation

$$\begin{aligned}\text{Percentage of moisture} &= \text{Wt. of moisture} / \text{Wt. of coal} \times 100 \\ &= (W_2 - W_3) / (W_2 - W_1) \times 100\end{aligned}$$

b. Determination of % Ash

Procedure:

1. First weigh the empty crucible (W_1 gm). Then take about 1gm of coal sample in the same crucible and weigh it again (W_2 gm).
2. Place the silica crucible in a muffle furnace at 750°C for 1 hour to complete the combustion.
3. Then the crucible is removed and allowed to cool in desiccators to room temperature

Downloaded from- www.abhijitgurav.wixsite.com/1234

For more information visit – www.youtube.com/c/chemistrylearners

and weigh it again (W_3 gm).

- The amount of residue remaining in the crucible corresponds to the ash content of the coal.

Observations

Observations	Value
Weight of empty crucible (W_1)	= gm
Weight of crucible with coal sample (W_2)	= gm
Weight of crucible with coal sample after heating (W_3)	= gm
Weight of moisture/ ash	$(W_3 - W_1) =$ gm

Formula and Calculation

Percentage of ash = Wt. of ash / wt. of coal sample $\times 100$

$$= (W_3 - W_1) / (W_2 - W_1) \times 100$$

c. Determination of % Ash

- Moisture free coal is heated upto 925°C temperature for 3-4 hrs.
- It is (crucible along with coal) then kept into desiccator for cooling.
- Again the weight of crucible is taken and % volatile matter is calculated as

Percentage of volatile matter = Wt. of volatile matter / wt. of dry coal sample $\times 100$

d. Determination of % Fixed carbon

This is actual carbon present in coal for combustion. It is calculated as

% of Fixed carbon = $100 - (\% \text{ Moisture} + \% \text{ Ash} + \% \text{ Volatile matter})$

Significance of proximate analysis

- Moisture lowers the calorific value and takes more the time to burn fuel. Therefore, less the % moisture, better is the quality of coal.
- Ash is non-reducible byproduct formed after burning of coal. It has no calorific value. Therefore, less the % ash, better is the quality of coal.
- Volatile matter elongates the flame size. Therefore, calorific value decreases. Therefore, less the % Volatile Matter, better is the quality of coal.

- d. Fixed carbon is giving calorific value. Therefore, more the % fixed carbon, better is the quality of coal.

Chemistry Learners